

How to use OR Info

1

Smart OR Info

Select one of the options below to continue

Select '**Check-In**' to add new instruments.

2

CHECK-IN

Select OR

Select an operation room and click '**Continue**'.

3

CHECK-IN • OR 1

Select Surgery

Choose the relevant surgery and click '**Start selected**'.

4

Hold each instrument near the RFID reader to register it.

5

Click the '**Disposables**' tab to view and add disposable items.

6

Click '**Save & Continue Later**' to save your check-in session.



Scan to Access
Quick Start Guide