


How to Pack a Case Cart

1



Scan Case Cart
Click to read the case cart for packing

Start Scanning

Enter Barcode Manually

Click **'Start Scanning'**.

2

Select Surgery for Case Cart #23

Search

In Progress for Current Case Cart

24/10/2024 12:30
Operation Room 1 - [446339] Pancreaticoduodenectomy with Portal Vein Reconstruction
2nd Floor, Main Hospital (Copenhagen) | Dr. Jens Peter

Upcoming surgeries

24/10/2024 10:45
Operation Room 2 - [446338] Laparoscopic Cholecystectomy
2nd Floor, Main Hospital (Copenhagen) | Dr. Jens Peter

Assigned

24/10/2024 09:25
Operation Room 3 - [446340] Coronary Artery Bypass Grafting (CABG) Surgery
2nd Floor, Main Hospital (Copenhagen) | Dr. Jens Peter

Cancel Continue

Select the surgery and click **'Continue'**.

3

Case Cart #23

Operation Room 1 (1st floor) | Operation Room 2 (2nd floor) View +1 more

Add manually Unlock

Search packed items Packed: 2/9

TYPE	SURGERY	STATUS	REMOVE
Tray	446339 Laparoscopic Cholecystectomy	Packed	
Tray	446339 Laparoscopic Cholecystectomy	Packed	
PE, 8 CH	446339 Laparoscopic Cholecystectomy	Pending 0/6	
Pouch	446333 Cataract Surgery	Pending	

To pack instruments, place each item on the RFID scanner.

4

Case Cart #23

Operation Room 1 (1st floor) | Operation Room 2 (2nd floor) View +1 more

Add manually Unlock

Search packed items Packed: 2/9

SURGERY	STATUS	REMOVE
446339 Laparoscopic Cholecystectomy	Packed	

If extra items / sets are needed, click **'Unlock'**. Otherwise, proceed to Step 7.

5

Unlock case cart

Unlocking the cart allows you to add items outside the recipe.

Reason for unlocking

Select

Notes

Add notes (optional)

Cancel Unlock

Select a reason for adding extra items, then click **'Unlock'**.

6

Add items manually

Select from List Enter Barcode

Search

Popular searches

- General Use Pouch 34155BF9C800194000004E2
- GUIDE NEEDLE, ANG., LANCET SHAPE, 8 CH KLS Martin - 200100
- GUIDE NEEDLE, ANG., LANCET SHAPE, 10 CH KLS Martin - 200101
- GUIDE NEEDLE, ANG., KNIFE SHAPE, 12 CH KLS Martin - 200102

Cancel Continue

Select additional disposables click **'Continue'**.

7

Case Cart #23

Operation Room 1 (1st floor) | Operation Room 2 (2nd floor) | View +1 more

Scanned Items (2) | Differences (7) | Add manually | Unlock | →

Scan barcode to add items | Search packed items | Packed: 2/9

FULL NAME	TYPE	SURGERY	STATUS	REMOVE
gyn_gen_uro_základ_maly A-01 34155BF9C800194000004E1	Tray	446339 Laparoscopic Cholecystectomy	Packed	×
gyn_gen_uro_základ_maly A-02 34155BF9C800194000004E2	Tray	446339 Laparoscopic Cholecystectomy	Packed	×

Scan | Config | Activity

Once all items are packed, click **'Finish'**.

8

Are you sure you want to finish packing Case Cart #23?

CURRENT SESSION Packing 20/20

Operation Room 1 (1st floor)
[446338] Laparoscopic Cholecystectomy

LAST PACKED 23/12/2024 13:01 Packing 20/20

Operation Room 2 (2nd floor)
[446337] Cataract Surgery

LAST PACKED 23/12/2024 12:45 Packing 20/20

Operation Room 3 (3rd floor)
[446339] Pancreaticoduodenectomy with Portal Vein Reconstruction

Cancel | Save & Continue Later | **Yes, Finish Packing**

Click **'Yes, Finish Packing'**.

9

Packing discrepancy detected (8/9)

Enter PIN code to force finish:

• • • • • •

Notes

Add additional notes (optional)

Cancel | **Force Finish**

In case of discrepancy, enter PIN and click **'Force Finish'**.

10

✓ Packing Completed

Case Cart #23 has been packed; please review the summary below:

Operation Room 1 (1st floor); Operation Room 2 (2nd floor);
Operation Room 3 (3rd floor)

gyn_gen_uro_základ_maly A-01	1/1	✓
gyn_gen_uro_základ_maly A-02	1/1	✓
GUIDE NEEDLE, ANG., KNIFE SHAPE, 8 CH	6/6	✓
General Use Pouch	1/0	✓

Summary	Sets	Pouches	Disposables	Total Items
	3/3	1/1	6/6	10

Print Label | View Report | **Back to Home Screen**

Packing is completed.



Scan to Access
Quick Start Guide